



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 3038



Lab Supplies, Disinfectant and Cleaner, 1 U.S. Quart, PREempt RTU

DESCRIPTION: This EPA-registered disinfectant utilizes AHP[®], a patented synergistic blend of commonly used ingredients that when combined with low levels of hydrogen peroxide dramatically increase its germicidal potency and cleaning performance. AHP[®] combines hydrogen peroxide with surfactants, wetting and chelating agents to create an ideal one-step disinfectant cleaner. Ideal for cleaning and disinfecting environmental surfaces in cleanrooms and laboratory areas including workstations, fume hoods, equipment and other hard non-porous environmental surfaces.

FEATURES:

- ✔ Contains 0.5% Hydrogen Peroxide
- ✔ 3 year shelf life (opened and unopened)
- ✔ Effective, fast-acting, one-step disinfectant cleaner
- ✔ Convenient, ready-to-use disinfectant for daily use
- ✔ Low residue and compatible with a wide range of materials
- ✔ Odor-free for good operator acceptability
- ✔ Emits no VOCs and does not require personal protection equipment (PPE)
- ✔ Hydrogen peroxide breaks down into water and oxygen for easy disposal and minimal, to no environmental impact
- ✔ Not classified as dangerous goods
- ✔ Contains 0.5% Hydrogen Peroxide
- ✔ 2 year shelf life

TECHNICAL SPECIFICATIONS:

Volume (U.S Quart/oz/L)	1 / 32 / 0.9
Dispensing Method	Manual pour
Scent	Odor-free
Concentration of Actives	0.5% Hydrogen Peroxide

DIRECTIONS OF USE:

Disinfection of non-porous hard surfaces (i.e. exterior of fume hoods, tabletops and other environmental surfaces):

1. Heavily soiled surfaces require cleaning prior to disinfection.
2. Apply to surface with cloth or disposable wipe.
3. Ensure surface remains wet for 3 minutes.

Disinfection of Non-Critical, Pre-cleaned Instruments and Medical Devices:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and dried.
2. Clean and rinse lumens of hollow instruments before filling with this product.
3. Using either tray or an ultrasonic unit, immerse instrument in PREempt™ RTU for 3 minutes.
4. Wipe with a clean, damp cloth or paper towel and allow to air dry.
5. Discard solution after each use.

EFFECTIVE AGAINST A BROAD-SPECTRUM OF VIRUSES, BACTERIA AND FUNGI IN 3 MINUTES!

Virucidal: 3 Min.

- ✔ Proven effectiveness against the Poliovirus Type 1, Sabin strain type 1 which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses
- ✔ HIV-1 Human Immunodeficiency Virus (HIV)
- ✔ Human Coronavirus 229E
- ✔ Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses
- ✔ This product has demonstrated effective against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

Bactericidal: 3 Min

- ✔ *Staphylococcus aureus*
- ✔ *Pseudomonas aeruginosa*
- ✔ *Salmonella enterica*
- ✔ Methicillin-resistant Staphylococcus Aureus (MRSA)
- ✔ Vancomycin-resistant Enterococcus faecalis (VRE)

Fungicidal: 3 Min.

- ✔ *Trichophyton mentagrophytes*

Tuberculocidal: 3 Min.

- ✔ *Mycobacterium terrae*

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Broad-Spectrum Sanitizing: 30 sec.

-  *Staphylococcus aureus*
-  *Pseudomonas aeruginosa*
-  *Salmonella enterica*
-  *Klebsiella pneumonia*
-  *Escherichia coli O157:H7*
-  *Campylobacter jejuni*

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM