



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6349**



## Adapter, Press-In Bottle, for 33 mm Bottle

**DESCRIPTION:** Press-in bottle adapters allow for increased dosing accuracy and fit securely in any bottle with the appropriate size opening. The Press-in Bottle Adapters allow a user to easily draw liquid with an oral dispenser, ensuring accurate dosing while avoiding spills. The adapters fit flush with the bottle opening so that the original cap can be placed on the bottle.

### Which Press-In adapter will fit my bottle?

<i>Find your bottle by closure code, outer thread diameter or inner neck diameter below.</i>			<u>Adapter Stock Code</u>
<u>Closure Code</u>	<u>Outer Thread Diameter</u>	<u>Inner Neck Diameter (ID)<sup>(2)</sup></u>	<u>Code</u>
20-xxx <sup>(1)</sup>	20 mm	15 mm	<b>8796</b>
24-xxx	24 mm	19 mm	<b>8795</b>
28-xxx	28 mm	22.5 mm	<b>8794</b>
33-xxx	33 mm	24 mm - 28 mm	<b>6349</b>

**Note:** MEDISCA cannot guarantee full standardization of all bottles on the market. Please test your adapter for final confirmation of fit.

(1) xxx refers to the thread finish, such as 400, 410 or 415.

(2) ID values are approximate.

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)