



# **QUICK JUST BECAME *FASTER!***



Now is the time to take your compounding practice to the next level by offering your patients MEDISCA's new Rapid Dissolve Tablet. Our innovative RDT Kit accommodates your customers' needs.

## > WHAT ARE RDTs?

Rapid Dissolve Tablets are designed to dissolve quickly on the tongue without any water making it an ideal alternative dosage form for patients who have difficulty swallowing. In short, they offer quick delivery of medication, rapid dissolution, and greater absorption in the blood.

## > WHAT'S INCLUDED IN THE RDT KIT




1

### MEDI-RDT BASE

The Medi-RDT Base is a **Sugar-Free**, finely granulated powder that is compatible with a wide range of active ingredients. Our unique RDT formula provides **significantly faster** wetting and disintegration times compared to similar products on the market:

- Wetting time: 5 sec. (versus 24 sec.)
- Disintegrating time: 93 sec. (versus 287 sec.)

Included in the kit: 100 g



2

### MEDI-RDT MOLD

Made out of an extremely durable, heat-resistant hospital grade polymer, the mold is designed for repeated and prolonged baking, up to 120°C (248F). Committed to providing the highest quality products and supplies, MEDISCA's Teflon-free molds were developed to avoid the risk of any particles potentially getting into the tablets and causing contamination.

The lateral containment system helps to avoid powder loss and facilitates the filling process, while the calibrated cavities allow you to compound tablets that comply with the average weight test.

Additional benefits include:

- Lightweight for easy handling and cleaning
- Dishwasher safe
- Design allows for a significant decrease in baking time

Included in the kit: 1 RDT mold with 96 cavities (750 mg)



*"The new RDT mold worked wonderfully! It is easier to fill, clean and it bakes for ½ the time of the others we have tried. In addition, the switch over was made without hassle, there was no need to change our current formula. We are very pleased with the results and excited that we have been able to reduce the compounding time."*

MATTHEW BERNSTEIN, R.PH., HM COMPOUNDING



## > WHY CHOOSE MEDISCA'S RDT KIT?

MEDISCA's Medi-RDT Sugar-Free Base has been put to the test against competitors and has proven to have superior performance. Furthermore, the Medi-RDT mold is made out of Teflon-Free heat resistant polymer which means that it is resistant to repeated and prolonged heating up to 120°C (248F). These are some of the reasons for which more and more compounding pharmacists are choosing the MEDISCA RDT Kit.



3

### MEDI-RDT BLISTER PACKS WITH SLEEVES FOR 160 DOSES

The customized air-tight blister packs include UV-resistant amber-colored sheets for optimal protection of light-sensitive drugs. Cold seal tamper evident foil laser labels for medication identification and sleeves are also included to be able to dispense to your patients.



4

### HOT HAND® PROTECTOR

With safety being a top priority, the Hot Hand® protector mitt provides superior protection when handling hot or cold objects. The gripping surface is studded with suction cups for a secure no-slip grasp to safely remove the Medi-RDT Mold from the convection oven.



5

### CUSTOMIZED FORMULA

MEDISCA NETWORK recognizes your need to get it right the first time which is why in order to ensure proper compounding method techniques for the Rapid Dissolve Tablet, MEDISCA NETWORK will provide you with a customized sample formula of your choice to help get you started with confidence.

6

### 50 MESH SIEVE AND RECEIVER PAN

Another essential tool included is the 50 mesh sieve and pan which serve to control particle size distribution as well as the flow property of the blend. This is to ensure content uniformity of the solid dosage form.





## ➤ GET STARTED TODAY!

Grow your practice and competitive edge with MEDISCA's comprehensive RDT Kit. We want to set you up for success as quickly as possible with all the necessary tools to elevate your compounding practice and surpass your patients' expectations with our Rapid Dissolve Tablet Kit.

MEDISCA wants you and your patients to experience how **QUICK JUST BECAME FASTER!**

### ADDITIONAL ITEMS SOLD SEPARATELY

- 1 MEDI-RDT MOLD (Product No.: 7311-01)**  
 Multi-dose mold with 3 cavity sizes – ideal for compounders requiring greater versatility and different dosage options.
  - 180 Cavities: 60x75 mg, 60x150 mg, 60x200 mg
- 2 MECHANICAL PREMIUM ISOTEMP CONVECTION OVEN (Product No.: 9011-01)**
  - 23.5"W x 25.75"D x 26"H, 70.8L / 2.5 cuft, 50-275°C / 122-527°F, Vertical Forced Air, w/One Shelf (Holds Up to 5)
- 3 CONVECTION OVEN (QUINCY) (Product No.: 8734-01)**
  - 20"W x 14.3"D x 25.5"H, 450°F / 232°C, Analog Temperature, Gravity Convection



CALL NOW OR VISIT OUR WEBSITE TO PLACE AN ORDER

🇨🇦 1.800.665.6334 | [WWW.MEDISCA.COM](http://WWW.MEDISCA.COM) | 1.800.932.1039 🇺🇸

**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

DISCLAIMER : For use in prescription pharmacy compounding only, in accordance with applicable law. MEDISCA makes no warranty or representation regarding the use of this product or fitness for any particular purpose, and any such use is entirely the responsibility of the pharmacist in the exercise of his/her professional expertise and judgment. Accordingly, in no event will MEDISCA be liable for any damages from the use or misuse of the product. Purchaser's sole remedy in the event of product defect or other damage, shall be the refund of the purchase price of the product.